# **Tests for Biological Molecules**

PREPARED BY: BESIR ZENELI

#### 1. Iodine Test For Starch

- •Summary: Purpose: Detecting the presence of starch.
- •Principle: Iodine reacts with starch, forming a blue-black color.
- •Observations: Brown to blue-black color change indicates the presence of starch.
- Application: Used in biology and food science.




## Animation – Iodine Starch

https://www.youtube.com/watch?v=YI-Imfp-f3g

#### 2. Biuret Test For Protein

- •Summary: Purpose: Identifying the presence of proteins.
- •Principle: Biuret reagent reacts with proteins, turning the solution purple.
- •Observations: Blue to purple color change indicates the presence of proteins.
- •Application: Important in biochemistry and food analysis.




## Animation — Biuret Protein

https://www.youtube.com/watch?v=pDfBOf3seUc

### 3. Benedict's Test For Glucose

- •Summary:Purpose: Detecting the presence of glucose.
- •Principle: Benedict's solution reacts with glucose when heated, changing color.
- •Observations: Blue to green, yellow, or brickred color change indicates presence of glucose.
- Application: Used in food industry and clinical chemistry.



## Animation – Benedict Glucose

https://www.youtube.com/watch?v=HxevnEdU9AI

### 4. Ethanol Emulsion Test For Fats

- •Summary:Purpose: Indicating the presence of fats.
- •Principle: Fats, being hydrophobic, form emulsions when mixed with ethanol.
- •Observations: Formation of a cloudy or milky emulsion.
- Application: Useful in food and lipid analysis.





Figure 9: The emulsion test for lipids.

## Animation – Emulsion Lipids

https://www.youtube.com/watch?v=wScbwdrlucA

#### 5. DCPIP Test for Vitamin C

- •Summary:Purpose: Detecting the presence of vitamin C (ascorbic acid).
- •Principle: DCPIP changes color when it reacts with vitamin C.
- •Observations: Color change from blue to colorless or pink.
- •Application: Applied in nutritional analysis and pharmaceuticals.



## Animation – DCPIP Vitamin C

https://www.youtube.com/watch?v=tnEYA8Nd4o8