



# Extracting DNA from a Banana

Prepared by: Besir Zeneli, BSc. in Biology

# Extracting DNA from a Banana at Home

A simple experiment using household materials

# Materials Needed

---

## List:

- 1/2 ripe banana
- Water
- Dish soap or liquid soap
- Cold alcohol (ethanol or rubbing alcohol)
- Sieve, coffee filter, or tissue
- Spoon or pestle
- Transparent glass or test tube
- Stick for collecting DNA



# Step 1 - Chop the Banana

---

Peel and finely chop **half a banana**.  
Place it into a cup for mashing.



shutterstock.com • 1349919032

## Step 2 - Mash the Banana

---

Use a spoon or pestle to mash the banana into a smooth paste.

Alternatively, you can use **plastic resealable freezer bag**.





## Step 3 - Prepare Lysis Solution

---

**Lysis** refers to the breaking down or destruction of cells, typically by the rupture of the cell membrane.

In another cup, mix 1 teaspoon of dish soap with 10 tablespoons of water.

This is called a lysis solution. Its function is to break the cell membrane of the cells.



## Step 4 - Mix the Lysis Solution

---

Pour the lysis solution into the mashed banana and stir well for at least 5 minutes.



## Step 5 - Filter the Mixture

---

Strain the mixture using a sieve, coffee filter, or tissue. Collect the liquid in another cup.





## Step 6 - Transfer to a Transparent Container

---

Pour the filtered liquid into a test tube or a tall, thin glass jar.



## Step 7 - Extract the DNA

---

Slowly pour cold alcohol along the side of the container. A white, stringy substance will appear – **this is the DNA**.



## Step 8 - Collect the DNA

---

Use a stick to lift the DNA strands out. Note that this is not pure DNA, as it contains proteins and other materials.



How to do it?

[https://www.youtube.com/watch?  
v=MYBuY5Nvg9s&t=67s](https://www.youtube.com/watch?v=MYBuY5Nvg9s&t=67s)